

Leipzig Applicator HDR End-to-End Testing using In-house Fabricated 16 mm EBT4 Radiochromic Film Discs

Nadia M. Vassell, MS, Charlotte Read, MS, Mario Serrano-Sosa, PhD, Sunshine Osterman, PhD, and Andrew Lukban, MS, Mount Sinai Health System, New York, NY

INTRODUCTION

Commissioning of horizontal Leipzig applicators for HDR brachytherapy typically includes determination of the effective distal source position and dose output verification. We previously reported an end-to-end testing methodology using optically stimulated luminescent dosimeters (OSLDs), in which measured doses were compared with published Monte Carlo reference data. Following the discontinuation of OSLDs, we developed a refined approach using in-house fabricated EBT4 radiochromic film discs. The purpose of this work was to evaluate whether EBT4 film can provide equivalent dosimetric verification for Leipzig applicator commissioning.

MATERIAL & METHODS

EBT4 radiochromic film discs (16 mm diameter) were fabricated in-house and calibrated in solid water using a 6 MV photon beam. Custom CAD-designed, 3D-printed holders were developed for both the film discs and horizontal Leipzig applicators to ensure reproducible alignment of the film along the applicator central axis. Expected doses to water at 3 mm depth were obtained from published Monte Carlo-based dose tables for horizontal Leipzig applicators with nominal diameters of 3 cm, 2 cm, and 1 cm. A flatbed scanner was used to create a dose calibration curve and measured doses were then compared with Monte Carlo reference values to evaluate dosimetric agreement.




Figure 1: Deconstructed CAD design showing 3D printed film holder, 16mm EBT4 film, alignment tool and Leipzig applicator set up.




Figure 2: EBT4 film calibration setup in solid water showing the 3D-printed alignment tool used to ensure reproducible positioning of the 16 mm film discs during irradiation.




Figure 3: Experimental set up for dose measurement of the 1cm horizontal Leipzig applicator, consisting of solid water, a 3D-printed Leipzig alignment tool and a film holder containing the 16 mm EBT4 film positioned at a depth of 3mm. Also shown are the 2 cm and 3 cm Leipzig applicators used in subsequent measurements.

RESULTS

Measured doses for all Leipzig applicator sizes were within 2% of published Monte Carlo reference values at 3 mm depth. Specifically, percent differences were 1.13%, 1.75%, and 0.94% for the 3 cm, 2 cm, and 1 cm applicators, respectively. These results also demonstrate agreement to that previously achieved using OSLDs.




Figure 4: EBT4 film calibration generated using 6MV irradiation. Data were fit using a 3rd order polynomial function with an R^2 of 0.9995 demonstrating excellent agreement between measured and fitted values.

Applicator/Size	Reading	Net Optical Density	Expected dose [cGy]	Measured dose [cGy]	Percent error [%]
Leipzig-3cm	71	0.065	193.45	191.25	1.13
Leipzig- 2cm	71	0.065	194.66	191.25	1.75
Leipzig-1cm	72	0.067	197.68	199.54	0.94

Table 1: Monte Carlo expected dose vs measured dose with 16 mm EBT4 film for horizontal Leipzig applicators with nominal diameters of 3 cm, 2 cm, and 1 cm.

SUMMARY / CONCLUSION

End-to-end testing of Leipzig applicators is challenging because their high-density construction limits the use of conventional CT-based, water-equivalent verification workflows within standard treatment planning systems. In-house fabricated 16 mm EBT4 radiochromic film discs demonstrated dosimetric agreement within 2% of published Monte Carlo reference values and provided an effective replacement for OSLDs. This methodology offers a reliable, reproducible, and clinically practical solution for end-to-end testing of Leipzig applicators in HDR brachytherapy commissioning.