

Improving the Accuracy and Detectability of HfO_2 Nanoparticles

Ramaswamy Sadagopan, MS, Sabrina Campelo, PhD, MS and Sandun Y. Jayarathna PhD.

Department of Radiation Physics

The University of Texas MD Anderson Cancer Center, Houston, Texas, USA

INTRODUCTION

High-atomic number hafnium oxide (HfO_2) nanoparticles have demonstrated potential as radio-enhancers, but their detectability using standard CT calibration is limited. This work aims to (1) establish CT calibration curves for HfO_2 solutions of known concentration and (2) evaluate dual-energy CT's capability to improve the minimum detectable concentration of standard CT.

METHODS

RMI (Gammex) phantom inserts held aqueous HfO_2 solutions, 0.1% to 10% by weight, imaged with a thoracic free-breathing protocol (120 kVp) and dual-energy (100/140 kVp Sn) with monoenergetic and kVp reconstructions. A water-filled-vial scan served as control.

Circular ROIs (6 mm radius) were centered in the upper vial liquid, above any settled material, and averaged over the central liquid slices. Solid-water regions gave the water-equivalent reference and noise.

RESULTS

In thoracic free-breathing images (120 kVp), mean HU increased monotonically with concentration, from 36.2 ± 25.5 HU at 0.1% to $2,750.5 \pm 138.8$ HU at 10% by weight.

Dual-energy monoenergetic images increased contrast per unit concentration; 50 keV VMIs read 33.6 ± 26.5 HU at 0.1% and $4,131.4 \pm 127.6$ HU at 10%. Across reconstructions, 0.5% was clearly detectable ($\text{CNR} > 5$) while 0.1% approached the noise floor ($\text{CNR} \approx 1$ to 2.5). Lower-energy VMIs gave the highest CNR.

CONCLUSIONS

- CT calibration curves for HfO_2 nanoparticles were successfully established and validated using a water-filled controls
- Dual-energy CT improved contrast per unit concentration compared with conventional 120 kVp CT.
- Low-energy virtual monoenergetic images provided the greatest contrast enhancement while maintaining image noise.
- HfO_2 concentrations $\geq 0.5\%$ by weight were clearly detectable, while 0.1% approached the detection limit.

These findings support the use of dual-energy CT for improved visualization of HfO_2 nanoparticles in radiotherapy applications.




Figure 1. CT calibration curve at 120 kVp (thoracic free-breathing) for HfO_2 nanoparticle solutions. Mean HU increases with concentration from 0.1% to 10% by weight. Error bars represent ± 1 SD (within-ROI noise).




Figure 2. Dual-energy CT calibration curves for HfO_2 solutions: virtual monoenergetic images at 50, 70, and 80 keV, plus 100 kVp and 140 kVp (Sn) reconstructions. Lower-energy VMIs give the greatest contrast per unit concentration. Error bars represent ± 1 SD.




Figure 3. Adjacent-concentration separability by reconstruction. Bars represent the separation index $z = |\Delta \text{HU}| / \sqrt{(\text{SD}_{\text{low}}^2 + \text{SD}_{\text{high}}^2)}$ between successive concentrations. Values above the dashed line ($z \geq 2$, $\sim 2\sigma$) are statistically resolvable, and the red "x" symbols flag reconstructions that fall below it for a given step.



Figure 4. RMI phantom with inserts and vials containing HfO_2 solutions of various concentrations.

INNOVATION

To the best of our knowledge, this is the first attempt to calibrate HfO_2 nanoparticle using clinical CT scanner.

**Presenting Author: rsadagop@mdanderson.org

*Corresponding Author: rsadagop@mdanderson.org