

Estimating Relative Levels of Linolenic Acid with Magnetic Resonance Spectroscopy of Fat Allylic Protons at 3 T

Iain Marchi¹ and Atiyah Yahya^{1,2}

1. Department of Oncology, University of Alberta, Edmonton, Alberta, Canada
2. Department of Medical Physics, Cross Cancer Institute, Edmonton, Alberta, Canada

ABSTRACT

Purpose: To exploit differences in J-coupling evolution of the ≈ 2 ppm fat allylic protons (neighbor the two ends of the double bond chain) to estimate relative linolenic acid content at 3 T. Oleic (mono-unsaturated), linoleic (di-unsaturated) and linolenic acid (tri-unsaturated) are the three primary unsaturated components of edible oils.

Methods: A 3 T Philips MRI scanner was used with a Point RESolved Spectroscopy (PRESS) sequence (first echo time, TE₁, 17 ms) to measure the response of fat allylic protons from oleic, linoleic and linolenic acid. TE values ranged between 40 ms and 300 ms in steps of 10 ms. A TE that resulted in relatively low yield for oleic and linoleic acid but a significantly higher output for linolenic acid was identified. Measurements were obtained with the identified TE and with a TE of 70 ms (resolves allylic resonance from that of the ≈ 2.2 ppm resonance) from oils of varying linolenic acid content (linseed, walnut, canola, soybean and peanut oil).

Results: A TE of 220 ms was found to be effective for identifying the presence of linolenic acid. J-coupling evolution resulted in inverted allylic signal for all three constituents. The allylic resonance areas obtained with a TE of 220 ms normalized to that acquired with a TE of 70 ms were -0.07, -0.22 and -2.8 for oleic, linoleic and linolenic acid, respectively. The higher signal yield for linolenic acid is likely a result of more J-coupling evolution rewinding due to weaker coupling. A correlation was found between linolenic acid contribution to total unsaturation and allylic signal area obtained with a TE of 220 ms normalized to that from a TE of 70 ms (R² = 0.99).


Conclusion: Allylic fat proton signal obtained with PRESS with a TE of 220 ms can provide an indirect means of estimating relative linolenic levels at 3 T.

INTRODUCTION


- Magnetic Resonance Spectroscopy (MRS) provides a non-invasive means of assessing fat composition¹.
- Fatty tissue in the body, like edible oils, consist of saturated and unsaturated fatty acids. The latter consists of oleic (mono-unsaturated), linoleic (di-unsaturated) and linolenic (tri-unsaturated) acid.
- In unsaturated fatty acids, there are two allylic proton groups (≈ 2 ppm) that neighbour the two ends of the double bond chain.
- In linolenic acid, the allylic group neighbours the methyl group (unlike in oleic and linoleic acid) resulting in differing J-coupling evolution².

Purpose: To demonstrate that the J-coupling evolution of the allylic protons can be exploited to provide an indirect means of assessing relative linolenic acid content at 3 T.

Oleic acid



Linolenic acid



METHODS AND MATERIALS

- A 3 T Philips MRI scanner was employed.
- A Point RESolved Spectroscopy (PRESS) sequence (first echo time, TE₁, 17 ms) was used to measure the response of fat allylic protons from oleic, linoleic and linolenic acid. TE values ranged between 40 ms and 300 ms in steps of 10 ms.
- A TE that resulted in relatively low yield for oleic and linoleic acid but a significantly higher output for linolenic acid was identified.
- Measurements were obtained with the identified TE and with a TE of 70 ms (resolves allylic resonance from that of the ≈ 2.2 ppm resonance) from oils of varying linolenic acid content (linseed, walnut, canola, soybean and peanut oil).

Oil	% saturated	% Oleic	% Linoleic	% Linolenic
Peanut	19.5	60.8	19.7	negligible
Canola	9.2	60.7	20	10.1
Soybean	16.6	21.7	53.4	8.3
Walnut	10.2	20.4	53.2	16.2
Linseed	10.7	24.4	8.3	56.6

RESULTS AND DISCUSSION

- Figure 1 shows the allylic resonance (right of the dashed line) acquired from oleic, linoleic and linolenic acid..
- A PRESS TE of 70 ms resolves the allylic resonance from that of the ≈ 2.2 ppm protons (left of the dashed line).
- A TE of 220 ms yields relatively low allylic signal area for oleic acid (-0.07 times that obtained with a TE of 70 ms) and linoleic acid (-0.22 times that with a TE of 70 ms) but yields significantly higher signal for linolenic acid (-2.8 times that with a TE of 70 ms).
- Figure 2 displays example spectra obtained with a TE of 220 ms from linseed, walnut and peanut oil (linoleic acid contribution to total unsaturation 63.4%, 18% and 0%, respectively).
- The signal yield is lower with decreasing linolenic acid content.
- Normalizing the allylic signal area to that obtained with a TE of 70 ms from each oil yields a linear correlation with linolenic acid content (Figure 3)
- Allylic fat proton signal obtained with PRESS with a TE of 220 ms can provide an indirect means of estimating relative linolenic levels at 3 T.

Figure 1




Figure 2





Figure 3



REFERENCES

- Ren *et al.*, *Journal of Lipid Research*, **49**, p. 255 (2008)
- N. Bly *et al.*, *European Journal of Lipid Science and Technology*, **127**:e202400081 (2025)

ACKNOWLEDGMENTS

The Natural Sciences and Engineering Research Council of Canada (NSERC) is gratefully acknowledged for providing both student (NSERC USRA) and supervisor grant support (NSERC Discovery).

Philips Healthcare is acknowledged for research support.

